



## Manual Microscope Stage EK 100 × 100 for Upright Microscopes

Specifications	
Travel range	max. 102 × 102 mm
Drive	coaxial drive / rack – pinion
Material	aluminium
Surface	anodic coating, black lacquered
Weight	approx. 1.6 kg
Load capacity	5.0 kg

### Order Information

Order No.: 00-22-401-0000

### Product Features

- ▶ for reflected and transmitted light applications
- ▶ exchangeable stage inserts for different kinds of applications

### Accessories

Glass plate (included in delivery)  
 Order No.: 00-24-401-0810

For more accessories please visit:  
[www.marzhauser.com](http://www.marzhauser.com)

