



▶ **Motorized Measuring Stage MT mot 50 x 50** for measuring in transmitted light

Measuring range:	50 x 50 mm (2" x 2")
Repeatability:	< 1 µm
Accuracy:	absolut accuracy in µm: measuring length / 40 +2,5
Resolution:	< 0,01 µm
Motor:	2-phase stepper motor
Max. travel speed:	60 mm/s*, with 1 mm ball screw pitch
Weight:	3,500 g
Load capacity:	5,000 g

Order Information

MT mot 50 x 50,
with 1 mm ball screw pitch
Order No.: 00-36-102-0000

Product Features

▶ for measuring in transmitted light

Accessory

Glass plate (included in delivery)
Order No.: 00-30-101-0810

For further accessories please refer to:
www.marzhauser.com

Dimensions:

