



Manual Micromanipulator MD 4 with second electrode holder

| Specifications | |
|-----------------|---|
| Travel range | X-axis: 37 mm Y-axis: 20 mm Z-axis: 25 mm |
| Fine adjustment | 10 mm |
| Material | aluminium |
| Surface | anodic coating, black lacquered |
| Weight | approx. 700 g |

Order Information

Order No.: 00-42-105-0000 (right-sided, with clamp 1/2")
 Order No.: 00-42-106-0000 (left-sided, with clamp 1/2")
 Order No.: 00-42-121-0000 (right-sided, with tilting device)
 Order No.: 00-42-122-0000 (left-sided, with tilting device)

Product Features

- ▶ manual movement in 3 axes (XYZ)
- ▶ second electrode holder is independently turn-/tiltable and shiftable in XY direction
- ▶ several micromanipulators can be set up side by side
- ▶ coarse adjustment with 0.1 mm reading accuracy
- ▶ fine adjustment with 0.01 mm reading accuracy (micrometer screw)
- ▶ optional Z-axis stop avoids lowering of the micromanipulator (drift) when working with weighty loads

Accessories

Z-axis stop
 Order No.: 00-45-201-6000

Magnetic stand M1 (turn-/tiltable, adjustable in height)
 Order No.: 00-49-901-6000

Magnetic stand M2 (turn-/tiltable)
 Order No.: 00-49-902-6000

