




1-Axis Joystick for TANGO Controllers

Specifications	
Number of axes	1 (Z)
Interface	HD D-Sub 15-pin male plug connection to TANGO HDI
Material/surface	control panel: membrane keyboard housing: ABS-plastic with haptic coating housing bottom: stainless steel
Dimensions	229 × 122 × 90 mm (L × W × H)
Weight	approx. 1.0 kg



Order Information

1-Axis Joystick (Z-axis)
 Order No.: 00-76-100-0823

1-Axis Joystick (alternative axis)
 Order No.: on request

Product Features

- ▶ sensitive, manual positioning
- ▶ various configuration options
- ▶ free and precise programming of travel speed
- ▶ control of the LED illumination
LED 100 possible
- ▶ 4 freely programmable function keys (e. g. for position capture, velocity select etc.)
- ▶ membrane keyboard with high durability
- ▶ soft touch and insensitive haptic coating on the side surfaces
- ▶ easy cleaning

For more accessories for microscopy please visit: www.marzhauser.com