



## 2-Axes Joystick for TANGO Controllers

### Specifications

|                  |   |
|------------------|---|
| Number of axes   | 2 (XY)  |
| Interface        | HD D-Sub 15-pin male plug connection to TANGO HDI   |
| Material/surface | control panel: membrane keyboard<br>housing: ABS-plastic with haptic coating<br>housing bottom: stainless steel |
| Dimensions       | 229 × 122 × 90 mm (L × W × H)   |
| Weight           | approx. 1.0 kg  |


### Order Information

2-Axes Joystick  
 Order No.: 00-76-200-0820

### Product Features

- ▶ sensitive, manual positioning
- ▶ various configuration options
- ▶ free and precise programming of travel speed
- ▶ control of the LED illumination  
LED 100 possible
- ▶ 4 freely programmable function keys (e. g. for position capture, velocity select etc.)
- ▶ membrane keyboard with high durability
- ▶ soft touch and insensitive haptic coating on the side surfaces
- ▶ easy cleaning

### Dimensions



For more accessories for microscopy please visit: [www.marzhauser.com](http://www.marzhauser.com)