




3-Axes Joystick Digital for TANGO Controllers

Specifications

Number of axes	3 (XYZ)
Interface	HDI digital (Mini USB)
Material/surface	control panel: membrane keyboard housing: ABS-plastic with haptic coating housing bottom: stainless steel
Cable length	2.9 m
Weight	approx. 1.0 kg

Dimensions



Order Information

3-axes joystick digital
 Order No.: 00-76-300-1820

Product Features

- ▶ sensitive, manual positioning
- ▶ various configuration options
- ▶ free and precise programming of max. travel speed
- ▶ control of the LED illumination
LED 100 possible
- ▶ 4 freely programmable function keys (e. g. for position capture, velocity select etc.)
- ▶ membrane keyboard with high durability
- ▶ soft touch and insensitive haptic coating on the side surfaces
- ▶ easy cleaning

Compatibility

The 3-axes joystick digital is compatible with all TANGO controllers with digital HDI interface (e. g. TANGO integrale, TANGO 3 mini).