



4-Axes Joystick with Multi-Function Wheel for TANGO Controllers

Specifications

Number of axes	4 (XYZA)
Interface	HD D-Sub 15-pin male plug connection to TANGO HDI
Material/surface	control panel: membrane keyboard housing: ABS-plastic with haptic coating housing bottom: stainless steel
Dimensions	229 × 139 × 96 mm (L × W × H)
Weight	approx. 1.0 kg


Order Information

4-Axes Joystick
 with Multi-Function Wheel
 Order No.: 00-76-400-0820

Product Features

- ▶ sensitive, manual positioning
- ▶ side-mounted multi-function wheel allows ergonomic operation
- ▶ various configuration options
- ▶ assignment of the multi-function wheel to any axis
- ▶ free and precise programming of travel speed
- ▶ control of the LED illumination
LED 100 with multi-function wheel possible
- ▶ 4 freely programmable function keys (e. g. for position capture, velocity select etc.)
- ▶ membrane keyboard with high durability
- ▶ soft touch and insensitive haptic coating on the side surfaces
- ▶ easy cleaning

Dimensions



For more accessories for microscopy please visit: www.marzhauser.com