




1-Axis Joystick for TANGO Controllers

Specifications	analog	digital
Number of axes	1 (Z)	
Interface	HD D-Sub 15-pin	Mini USB
Compatibility	TANGO Desktop TANGO PCI-E	TANGO Desktop neo TANGO PCI-E neo TANGO 3 mini ¹ TANGO integrale ¹
Material/surface	control panel: membrane keyboard housing: ABS-plastic with haptic coating housing bottom: stainless steel	
Dimensions	229 × 122 × 90 mm (L × W × H)	
Weight	approx. 1.0 kg	

Dimensions



Order Information

1-Axis Joystick analog
 Order No.: 00-7604-0100-0000

1-Axis Joystick digital
 Order No.: 00-7604-1100-0000

Joystick variants for alternative
 single axes (X, Y, A) available on request

Product Features

- ▶ sensitive, manual positioning
- ▶ various configuration options
- ▶ free and precise programming of maximum travel speed
- ▶ control of the LED illumination
LED 100 possible
- ▶ 4 freely programmable function keys
(e. g. for position capture, velocity select etc.)
- ▶ membrane keyboard with high durability
- ▶ soft touch and insensitive haptic coating on the side surfaces
- ▶ easy cleaning

For more accessories for microscopy
 please visit: www.marzhauser.com