


Analog Micrometer, illustration exemplary
 Order No. 00-30-000-6804




Digital Micrometer, Order No. 00-99-901-4889

Micrometer Screws for Manual Measuring Stages

Dimensions



Analog Micrometer, illustration exemplary
 Order No. 00-30-000-6804



Digital Micrometer, Order No. 00-99-901-4889

Product Features

- ▶ Analog
 - ▶ Lead screw pitch: 0.5 mm
 - ▶ Accuracy according to DIN 863-1
 - ▶ Threads of ground hardened steel
- ▶ Digital
 - ▶ Digital step: 0.00005 inch / 0.001 mm
 - ▶ Max. permissible error: ± 0.00015 inch / ± 0.0038 mm
 - ▶ Non-rotating spindle

For more products please visit:
www.marzhauser.com

Measuring range	Display	Resolution	Thimble	Mounting shaft	Order No.
25 mm	analog	10 μ m	20 mm	12 mm	00-30-000-6801
25 mm	analog	5 μ m	30 mm	12 mm	00-30-000-6802
50 mm	analog	10 μ m	20 mm	12 mm	00-30-000-6804
50 mm	analog	5 μ m	30 mm	12 mm	00-30-000-6805
50 mm	digital (inch/metric)	1 μ m	–	18 mm	00-99-901-4889