Märzhäuser Wetzlar GmbH & Co. KG
In der Murch 15 | 35579 Wetzlar | Germany
Tel.: +49 6441 9116-0 | Fax: +49 6441 9116-40
info@marzhauser.com | www.marzhauser.com





Fig. shows
Order No. 00-0023-0827-0000

Manual Microscope Stage EK 200 × 200 for Upright Microscopes

Specifications

Travel range	205 × 205 mm (8" × 8")
Drive	self-locking coaxial drive, friction drive with quick adjustment
Material	aluminium
Surface	anodic coating, black lacquered
Weight	00-0023-0817-0000 approx. 7.3 kg 00-0023-0827-0000 approx. 8.7 kg 48-0023-0807-0000 approx. 8.1 kg
Load capacity	max. 5 kg

Dimensions



Supporting legs not shown;
guide handle length varies depending on version

Order Information

Supporting legs	Order No.
straight, 90 mm	00-0023-0817-0000
L-shaped, 130 mm	00-0023-0827-0000
L-shaped, 90 mm	48-0023-0807-0000

Product Features

- precise fine positioning via self-locking coaxial drive
- quick adjustment allows for effortless sample positioning
- straight or L-shaped supporting legs

Accessories

Metal plate 246 × 246 × 6 mm
(included in delivery)
Order No.: 00-24-801-0011

Glass plate 246 × 246 × 4 mm
Order No.: 00-24-801-0810

For more products please visit:
www.marzhauser.com

TÜVNORD
TÜV NORD CERT GmbH
DIN EN ISO 9001
tuev-nord.de

Specifications are subject to change without notice.
Issue: October 2026